

Stratedigm 4-laser configuration

Filters	405 nm	488 nm	552 nm	640 nm
445/60	PacBlue\DAPI			
530/30	Aqua\BV510	FITC\Alexa488		
580/30	PacOrange	PE		
615/30	BV605*	PE-TxR\PE-dazzle	mCherry	
676/29	BV650	PerCP		APC\Alexa647
710/40	BV711	PE-Cy5.5		Alexa700
740 LP	BV786	PE-Cy7		APC-Cy7

*Note: violet and green lasers are collinear, so only a violet or green excited fluorochrome can be used in combination on the same detector. (ie. Can't do BV605 with mCherry)