

Stratedigm 3-laser configuration

Filters	405 nm	488 nm	640 nm
445/60	PacBlue\BV421		
530/30	Aqua\BV510	Fitc\Alexa488\GFP	
580/30	PacOrange	PE	
615/30	BV605	PE-TxR\PE-Dazzle	
676/29	BV650	PerCP	APC\Alexa 647
690/40	BV711	PerCP-Cy5.5/PE-Cy5.5	
740 LP	BV786	PE-Cy7	APC-Cy7

Note: Red and violet lasers are collinear, so only a red or violet excited fluorochrome can be used in combination with a blue excited fluorochrome on the same detector.